
DESCRIPTION OF PARTS

DATA

Min. measuring range: 35,6 °C - 38,4 °C
 Resolution: 0,1 °C
 Accuracy: +/- 0,1 °C
 Operating conditions: Temperature: -15 °C to 39 °C

⚠ It is only possible to measure temperatures that are above the ambient temperature!

Storage conditions: Temperature: -15 °C to 39 °C
 Purpose of use: Analog ovulation / basal thermometer for cycle control
 This medical liquid glass thermometer is suitable for determining the basal temperature for natural family planning in the private sector. It measures the smallest changes in the temperature profile of the female cycle to determine the best days of conception.
 Analog maximum thermometer with Galinstan-filling
 The manufacture and testing of the thermometers is conducted in compliance with the standards

- of EU Council Directive 93/42/EEC concerning medical devices
- of EN 12470-1 – Clinical thermometers, Part 1: Metallic liquid-in-glass thermometers with maximum device

Geratherm® is certified in accordance with Council Directive 93/42/EEC, and DIN EN ISO 13485 and is entitled to affix the CE-mark CE 0197 (Notified Body: TÜV Rheinland LGA Products GmbH).

WHAT IS THE BASAL TEMPERATURE?

Women who want to gain knowledge about their fertility should apply the basal (morning) temperature method. Subtle variations in the basal provide information on the most fertile days of the woman.

⚠ PRECAUTIONARY MEASURES

1. The Geratherm basal is expressly not a contraceptive.
2. It is only possible to measure temperatures that are above the ambient temperature!
3. Do not drop the thermometer or subject it to sudden impacts.
4. Do not bend or bite the tip of the thermometer.
5. Keep the thermometer out of children's reach.
6. Thermometers with a damaged outer casing must not be used (risk of injury).
7. The fluid galinstan must not be brought into contact with light metals or precious metals because galinstan reacts with these metals.
8. Please be careful while measuring to avoid injury.
9. The use of the thermometer is no substitute of a doctor's visit.
10. Always clean and disinfect the thermometer before and after use.

INSTRUCTION FOR USE

1. Before taking the patient's temperature, check whether the column of liquid is below the first numbered line (36 °C).
2. Clean the glass thermometer before and after each use. In doing so, please pay attention to the "Notes on Cleaning" below.
3. Position the thermometer on the body according to the method of measurement to be used (axillary, rectal or oral).
4. After about 4 min the thermometer will have reached the body temperature. The temperature can be read off by rotating the thermometer gently to and fro. The measuring fluid in the glass capillary of the thermometer shows the measured body temperature on the temperature scale.
5. After measuring, the liquid must be returned to the bottom of the scale. To do this, hold the thermometer as close as possible to the upper end and shake it with 15 to 20 rapid wrist movements.

WHEN AND WHERE SHOULD I CHECK MY TEMPERATURE?

- Determine your temperature every morning at the same time before you get up. Don't eat before measuring. Avoid physical exertion. With an insufficient sleep break (at least 5 to 6 hours) do not enter a measured value, skip the measurement.
- Always measure on the same measuring point within a cycle: orally, vaginally or rectally. Do not change!
- Enter the temperature value in the included table.
- For problems or menstrual disorders, you should consult your doctor.

VAGINAL - measurement time: min. 4 minutes
 Please check with your gynecologist about this method of measurement.

⚠ STOP if you meet any resistance.

ORAL USE - measurement time: min. 4 minutes

Doctors prefer oral measurement. It provides fast and accurate results. Place the thermometer bulb (probe) in one of the two heat pockets to the left or right side of the tongue.

RECTAL USE - measurement time: min. 4 minutes

Rectal measurement is preferred in the case of children. Insert the thermometer bulb (probe) into the anus (approx. 2 cm). You may use a little skin cream or baby oil on the end of the vessel.
 You may use a probe cover for the thermometer. Insert thermometer into the probe cover. Lubricate the probe cover with a water soluble jelly for easier insertion. Do not use petroleum jelly. Insert tip of probe no more than 2 cm into rectum.
⚠ STOP if you meet any resistance.

AXILLARY USE - measurement time: min. 4 minutes, normal temperature: approx. 36,7 °C

The axillary measurement method provides less accurate measurements than oral and rectal measurement.
 Wipe armpit with a dry towel. Place probe in the patient's armpit and keep the patient's arm pressed firmly against their side.

HOW DO I RECORD MY TEMPERATURE AND WHY IS THAT IMPORTANT?

The days of your cycle 1 to 42 are plotted on top of the table, you can add more if necessary. Above each day of the cycle you write each date. Perpendicularly, the temperatures of 35 ° C to 38 ° C are plotted. Draft your temperature by means of a point mark. Mark the days on which you have sexual intercourse.
 1-2 days after ovulation (ovulation), the basal body temperature increases by 0.2 to 0.5 ° C and remains at this high position until the beginning of the next cycle. (Outside the ovulation is common to observe a decrease by about 0.1 ° C.) In the event that ovulation has occurred, the temperature rise should be entered within 48 hours. When you have filled about 3 tables, you will receive a Ovulations scheme and recognize your individual fertility scheme.

HOW TO BEGIN?

The first day of your menstrual flow is day 1 of your BBT chart. Your temperature should drop when your menstrual flow starts, do not include spotting prior to your period as your first day. Record your temperature throughout your period. Make certain you note the actual day of the month in the space provided on your chart. Note any premenstrual symptoms and any illnesses or fever and start a new table when you get your period again.

WHAT IS OVULATION?

Ovulation occurs when the egg is released from the ovary. Certain changes in the body may occur, such as the presence of thin clear vaginal mucus that comes from the cervix just before ovulation, around the 10th to 12th day of your cycle. This fluid allows sperm to remain alive for days until ovulation. Other symptoms may include a slight pain on the right or left side at the time of actual ovulation. The most likely days to become pregnant, with a 28-day cycle are the 10th and 14th day. Knowing when you are ovulating can be helpful in determining when to avoid sexual intercourse or increase the likelihood of becoming pregnant.

GALINSTAN – A NEW NON-TOXIC MEASURING FLUID

The mercury-free thermometer, Geratherm basal, has been developed by Geratherm Medical AG and contains a patented liquid alloy consisting of gallium, indium and tin – called Galinstan.
 Galinstan is completely non-toxic and environment-friendly.

PHYSICIAN CONSULTATIONS

It is strongly recommended that you consult with your personal physician when beginning family planning procedures. They can further assist you with interpreting the information on your tracking chart.

NOTES ON CLEANING

The thermometer has a hermetically sealed, hygienic glass casing and can be disinfected without restriction.

- The temperature of the cleaning and disinfectant solution must not exceed 39 °C. When measuring, please take the necessary care to avoid any danger of injury.
- The cleaning of the glass thermometer before and after use can be done with cold soapy water. Disinfect the thermometer using a piece of cotton wool or cloth that has been soaked in alcohol, or dip the thermometer in a disinfectant solution.
- If the thermometer breaks, any spilled liquid can be removed using a brush, tissue or cloth dampened in soapy water. The broken pieces may be disposed of as household waste.


LIMITED WARRANTY

Subject to normal use, this thermometer is unlimited guaranteed against function and accuracy (except breakage of glass) on the part of the manufacturer from the date of purchase. All parts of the thermometer are covered by this warranty. Damage to the thermometer caused by improper treatment is not covered by the warranty.
 The thermometer is valid indefinitely – calibration is not necessary!

SYMBOL INDEX

	Follow the instructions		Caution		Manufacturer
	Fragile: handle with care		Store in a dry place		Quantity
	Max. and min. temperatures: Store between -15 °C and +39 °C		Batch code (mm/yyyy, month/year) on the reverse of the outer casing (applied by means of laser without LOT symbol)		

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