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**GB Geratherm® pH-balance**  
 pH test for the testing of natural female intimate acidity protection for self-testing

**USAGE INSTRUCTIONS** Please read carefully.

**PACKAGE CONTENTS**

- 5 hygienically packaged test applicators for the measurement of vaginal pH
- 1 user manual
- 1 pH color scale

**FOR YOUR INFORMATION**

An acidic vaginal pH value of 3.8 to 4.5 is a basic requirement for the optimal functioning of the body's own system of protecting the vagina. This system can effectively avoid colonization by pathogenic germs and the occurrence of vaginal infections. The most important and most natural protection against vaginal problems is therefore a healthy vaginal flora. There are many factors that regulate the composition and dynamics of normal vaginal microbial ecosystem and that maintain the natural protective functions of the vagina. One important factor is lacto bacteria, see diagram 1.

Under the influence of lacto bacteria, the glycogen released from the upper vaginal epithelium degrades to form the sugar glucose and then lactic acid. This results in an acid environment forming in the vagina, with a pH value that is usually 3.8 to 4.5. This acidic vaginal environment forms effective protection against infections, because the vast majority of bacteria cannot survive under such conditions. The physiological vaginal pH level is therefore dependent on the quantitative presence of the lactic acid.

The pH level in vagina is subject to fluctuations. Possible causes of an alteration in the vaginal pH level are:

- Bacterial vaginosis (typical bacterial colonization of the vagina)
- Bacterial mixed infections
- Sexually transmitted diseases
- Premature rupture of fetal membranes
- Estrogen deficiency
- Postoperative infected wounds
- Excessive intimate care

Diagram 1: Diagram showing the degradation of glycogen into lactic acid

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**WASH YOUR HANDS BEFORE STARTING THE TEST.**

**Open the package, remove a test and open the foil. Only touch the applicator in the round handle area.**

**Hold the applicator in the handle area between your thumb and index finger, as shown in diagram 1.**

**The thumb should always be kept on the side of the vaginal applicators opposite to the pH measurement zone. If you use the vaginal applicator in this way, each test will be performed correctly.**

**Avoid the pH measurement zone coming into contact with other liquids.**

**Sit with your knees spread or with one foot on the toilet seat or other firm surface.**

**Use your your free hand to gently spread the vaginal opening, carefully separating the outer labia.**

**Insert the vaginal applicator gently into the vagina, in the same way that you would use a tampon, see diagram 2, (approx. 1 to 3 cm into the vagina - but not further than the round handle area).**

**Inserting the vaginal applicators is completely pain-free.**

**Now press the instrument gently against the back wall of the vagina (with very light pressure) for approx. 10 seconds (pushing towards the rectum), so that the pH measurement zone becomes well moistened. Sufficient moisture is very important for the test evaluation.**

**Now carefully remove the vaginal applicator from the vagina.**

**READING AND EVALUATING THE TEST RESULTS**

Read the results of the extraction of the vaginal applicators, while the measurement zone is still damp. Compare the color in the pH measurement zone with the color scale and note the value on the color scale that is most similar to the pH level you have measured. After evaluation, dispose of the used applicator.

**Wash your hands before starting the test.**

**Open the package, remove a test and open the foil. Only touch the applicator in the round handle area.**

**Hold the applicator in the handle area between your thumb and index finger, as shown in diagram 1.**

**The thumb should always be kept on the side of the vaginal applicators opposite to the pH measurement zone. If you use the vaginal applicator in this way, each test will be performed correctly.**

**Avoid the pH measurement zone coming into contact with other liquids.**

**Sit with your knees spread or with one foot on the toilet seat or other firm surface.**

**Use your your free hand to gently spread the vaginal opening, carefully separating the outer labia.**

**Insert the vaginal applicator gently into the vagina, in the same way that you would use a tampon, see diagram 2, (approx. 1 to 3 cm into the vagina - but not further than the round handle area).**

**Inserting the vaginal applicators is completely pain-free.**

**Now press the instrument gently against the back wall of the vagina (with very light pressure) for approx. 10 seconds (pushing towards the rectum), so that the pH measurement zone becomes well moistened. Sufficient moisture is very important for the test evaluation.**

**Now carefully remove the vaginal applicator from the vagina.**

**READING AND EVALUATING THE TEST RESULTS**

Read the results of the extraction of the vaginal applicators, while the measurement zone is still damp. Compare the color in the pH measurement zone with the color scale and note the value on the color scale that is most similar to the pH level you have measured. After evaluation, dispose of the used applicator.

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