
<b>Adizivisches Digital-Ag Fieberthermometer Geratherm® solar speed GT-161/1</b>	<b>Geratherm Medical-Ag Deutschland</b> CE0197	<b>(D)</b>
----------------------------------------------------------------------------------	---------------------------------------------------	------------

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
## Instruction for use

Please read carefully before using.

SPECIFICATIONS	
Range:	32.0 ~ 43.9 °C Below 32.0 °C, display L °C Above 43.9 °C, display H °C
Intended use:	Measurement of body temperature
Resolution:	0.1 °C
Accuracy:	±0.1 °C (34.0 °C to 42.0 °C) ±0.2 °C (other ranges)
Display:	Liquid crystal display 3 1/2 digits
Energy supply:	Integrated solar cell
Power consumption:	0.15 milliwatts in measurement mode
Power capacity:	Approximately 25 temperatures per 60-minute recharge
IP20 protection against foreign bodies < 12,5 mm)	
Dimensions:	147 mm x 21.3 mm x 15 mm
Weight:	Approx. 18.2 grams
Beeper:	Approx. 8 secs. sound signal when peak temperature reached
Memory:	Built-in memory enabling display of the previously measured value
Shut-off:	Automatic shut-off in abt. 3 minutes
Operating Conditions:	Temperature: -10 °C ~ 40 °C Humidity: 15 % ~ 95 % R.H.
Storage Conditions:	Temperature: -25 °C ~ 55 °C Humidity: 15 % ~ 95 % R.H.
Guarantee of Quality:	EN 12470-3 Clinical thermometers - Part3: Performance of compact electrical thermometers with maximum device. Geratherm® is certified in accordance with Council Directive 93/42/EWG and EN ISO 13485 and is entitled to affix the CE-mark CE0197 (Notified Body: TÜV Rheinland LGA Products GmbH).

Electromagnetic compatibility (EMC):
The present medical product is a device with extremely high degree of noise immunity. Therefore, the enclosure with the EMC tables was not printed out.

## ! PRECAUTIONS

- Use of this instrument is not intended as a substitute for consultation with your physician. Prolonged high fever requires medical attention especially for your children. Please contact your physician.
- Child's temperature shall be taken under adult's surveillance.
- Keep out of the reach of unsupervised children. Choking resulting from a child swallowing a small part that has become detached from the thermometer.
- Do not boil, bite, bend, drop or disassemble this thermometer.
- Do not attempt to disassemble the unit.
- Avoid taking temperature within 30 minutes after physical exercises, bathing or dining.
- Do not store, use and transport the unit under direct sunlight, at a high temperature, or in high humidity or dust.
- Performance may be degraded.
- Stop using the thermometer if it operates erratically or if the display malfunctions.
- If the thermometer has been stored at below-freezing temperatures, allow the thermometer to warm naturally to room temperature before using.

POWER RECHARGE
1. Before being used for the first time, the thermometer must be charged in good light for at least 5 hours.
2. When the low power symbol, <span>☐</span> , appears, the solar cell needs to be recharged at least 60 minutes.
3. Turn the power off before recharging.
4. The solar cell should face the light source during recharging.

## HOW TO USE

- Clean or disinfect the probe prior to each measurement. Position the probe orally, rectally or axillary (see measurement methods).
- Press the ON/OFF button. The display will read **1888**.
- Release the ON/OFF button and the display will show L °C with °C flashing.
- Place the probe in the appropriate position (oral, axillary or rectal).
- As soon as the thermometer does not register any further relevant increase in temperature, °C stops flashing and the maximum value measured up to that time is shown, accompanied by a series of audible signals.
- The thermometer continues measuring the temperature for a further minute after the audible signals begin, and in this time there may possibly be a further slight increase in the temperature shown. After that, the maximum temperature measured is recorded and stored; it remains visible on the LCD display for 3 minutes until the thermometer is automatically switched off. Or short press the ON/OFF button to switch off directly.
- The unit will automatically turn off in 3 minutes (approx.). However, to prolong battery life, it is best to turn the thermometer off by pressing the ON/OFF button once the temperature has been noted.

## ! NOTE:

You may put a sanitary plastic probe shield on the sensor tip to help prevent cross-infection and recontamination. Use of such probe cover may result in a 0.1 °C difference from actual temperature. Dispose of any used probe shield and store the unit in a safe place.

## SPECIAL FEATURES

**Normal Temperature Alarm:**
An alarm will sound when the maximum temperature has registered and is ready-to-read. The „Normal Temperature“ alarm (four quick beeps followed by a pause) will sound repeatedly for approximately 8 seconds.
**Fever Alarm:**
If the temperature registers higher than 37.5 °C, the fever alarm (a quick beeping sound which occurs continuously for about 8 seconds) will sound when the peak temperature has registered.
**Last Memory Recall:**
When the ON/OFF button is depressed to activate the thermometer, display of 188.8 will change to show the last measured temperature then return to „ready-to-measure“ mode and the recalled temperature will be erased.

### MEASUREMENT METHODS

! It is only is possible to measure temperatures that are above the ambient temperature.

### ORAL USE (about 45 - 90 seconds)

Place the probe well under the patient's tongue. Instruct the patient to keep their mouth closed.

A normal temperature by this method is usually considered between 36. 0 °C and 37.5 °C.

### RECTAL USE (about 30 - 60 seconds)

If rectal temperature is recommended by your doctor, you may use a probe cover. Insert thermometer into the probe cover. Lubricate the probe cover with a water soluble jelly for easier insertion. Do not use petroleum jelly. Insert tip of probe no more than 1.5 cm into rectum. STOP if you feel resistance. In general, the rectal temperature is around 0.5 °C to 1.0 °C **higher** than oral temperature.

### AXILLARY USE

Wipe armpit with a dry towel. Place probe in the patient's armpit and keep the patient's arm pressed firmly against their side. In general, the axillary temperature is around 0.5 °C to 1.0 °C **lower** than oral temperature. Note that in order to achieve better axillary temperature measurement result, **a minimum measuring time of 2 to 3 minutes is recommended regardless of the beep sound**.

! If you compare temperatures between two thermometers, it is possible to have two different readings even if both thermometers are within their accuracy tolerances. Differences in how you take each reading and the thermometer's tolerance range can account for nearly 0.1 °C variability.

**REMEMBER: Find your normal temperature and watch for changes.**

## CLEANING AND DISINFECTION

For household use, always clean the probe head before and after each use.

Wash the tip with a solution of mild detergent and cool water. Disinfect the thermometer by wiping the sensor and all parts are covered by this warranty excluding the battery.

The warranty does not cover damage to your unit due to improper handling. We recommend that the accuracy of the unit be checked every 2 years by an authorized laboratory.

## DISPOSAL

Observe the applicable regulations when disposing of the thermometer and batteries. This thermometer must not be disposed of together with domestic waste. All users are obliged to hand in all electrical or electronic devices, regardless of whether or not they contain toxic substances, at a municipal or commercial collection point so that they can be disposed of in an environmentally acceptable manner. Please remove the batteries before disposing of the thermometer. Do not dispose of old batteries with your household waste, but at a battery collection station at a recycling site or in a shop.

## LIMITED WARRANTY

This thermometer is guaranteed for two years from the date of purchase against manufacturer's defect under normal use. If your unit does not function properly due to defective parts or assembly we will repair it free of charge.

All parts are covered by this warranty excluding the battery. The warranty does not cover damage to your unit due to improper handling. We recommend that the accuracy of the unit be checked every 2 years by an authorized laboratory.

SYMBOL INDEX					
	Follow instructions for use		Type BF applied part		Manufacturer
	Keep dry		Caution! Read the instruction manual		Batch code (mm/yyyy, month/year)
	Store between 15 <span> </span> % and 95 <span> </span> % R.H.		Store between -25 <span> </span> °C and +55 <span> </span> °C		Geratherm Medical AG Fahrenheitsstrasse 1 98716 Geschwenda Germany www.geratherm.com
	The device must not be disposed of with household waste.		Numero del lotto (mm/aaaa; mese/anno)		Geratherm Medical AG Fahrenheitsstrasse 1 98716 Geschwenda Germany www.geratherm.com

<b>Medizinisches Digital-Ag Fieberthermometer Geratherm® solar speed GT-161/1</b>	<b>Geratherm Medical-Ag Allemagne</b> CE0197	<b>(F)</b>
-----------------------------------------------------------------------------------	-------------------------------------------------	------------

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