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CE0197


**Instructions**

*Self-testing kit for the diagnosis of urinary tract infections*

**SPECIFICATIONS**

The Geratherm® URINARY TRACT INFECTION TEST consists of a test strip with 4 reaction fields. These can quickly check for the following parameters: leukocytes, nitrite, blood and protein in the urine.

**REAGENTS AND DETECTION**

Reagent	Test duration	Composition	Characteristics
Leukocytes (LEU)	2 minutes	Amino acid ester pyrazole derivative, diazonium salt, buffer solution, non-reactive ingredients.	Detects leukocytes from Levels of 1 - 1w5 white corpuscles/Leu/mL in urine.
Blood (BLO)	1 minute	3,3',5,5'-tetramethylbenzidine (TMB), di-isopropylbenzene dihydroperoxide, buffer solution, non-reactive ingredients.	Detects haemoglobin from levels of 0.018 - 0.060 mg/dl or 5 - 10 Ery/µl in urine samples with an ascorbic acid content of <50 mg/dl.
Nitrite (NIT)	1 minute	P-Arsanilic acid, N-1(1-naphthyl)-ethylene-diamine, non-reactive ingredients.	Detects sodium nitrite from levels of 0.05 - 0.1 mg/dl in urine with low density and less than 30 mg/dl ascorbic acid.
Protein (PRO)	1 minute	Tetrabromphenol blue, buffer solution, non-reactive ingredients.	Detects protein from 7.5 - 15 mg/dl (0.75 - 0.5 g/L).

**PRECAUTIONS**

**Please read the instructions carefully before conducting the test.**

- For home use.
- Only to be used as a urine test, not as a blood test.
- Do not use after the expiry date.
- Keep out of children's reach.
- For in-vitro diagnosis, not for internal application.
- After use, the test kit must be disposed of in compliance with local waste disposal regulations.

**STORAGE AND USE**

- Store in a dry place at a temperature between 2 ° and 30 °C.
- Do not freeze; do not expose to direct sunlight.
- Keep the test strips in the sealed pouch until ready to use.
- If the test strips are exposed to the air for a long period, the reagents may be destroyed, leading to false test results.
- Do not use the test strips if the pouch is torn or damaged.
- Please do not touch the reaction fields on the test strips.
- If the test strips are faded, discard them; they may be defective.


**SCOPE OF DELIVERY**

Materials provided		
• Urine test strips	• Colour chart	• Instructions
Materials required but not provided		
• Sample collection container	• Clock or timer	• Kitchen paper

**INSTRUCTIONS FOR USE**

- Before testing ensure that the test strip and the urine sample are at room temperature (15 - 30 °C). It is best to use the first urine of the morning and to take the sample midstream. Please take the necessary hygiene measures before conducting the test.**
- Do not remove the test strip from the aluminium pouch until just before use, and grasp it as far away from the reaction fields as possible.
  - Do not hold the test strip in the urine stream until it has been flowing for 1-2 seconds. Ensure that all the reaction fields are thoroughly immersed in the urine (see Fig. 1).

- Note:**  
**You can also collect the urine in a clean, dry container.**  
In this case, dip the reaction fields completely in the urine and remove the test strip again immediately in order to prevent the reagents dissolving.  
When removing the strip from the container, carefully wipe it on the edge of the container to remove excess urine (see Fig. 1).

- Hold the test strip horizontal and lay it on a piece of absorbent tissue (e.g. a paper handkerchief) in order to prevent the chemical substances on the reaction fields and/or dirty hands coming into contact with the urine (see Fig. 2).
- Carefully compare the reaction fields with the corresponding colour bars on the colour chart. To do this, hold the test strip next to the colour bars (see Fig. 3).

**TEST PRINCIPLE AND EXPECTED VALUES**

**Leukocytes:**  
This test can be used to detect the presence of granulocyte esterases. The esterases break down a derivatised pyrazole amino acid ester, thus releasing hydroxy-pyrazole. This pyrazole reacts with a diazonium salt and produces colours ranging from pinkish beige to violet. Normal urine samples generally provide negative test results. If traces of leukocytes are present, this raises questions as to their clinical significance. If leukocytes are detected, we recommend that the test be repeated using a fresh urine sample from the same patient.

**Blood:**  
This is based on the pseudo-peroxidative activity of haemoglobin, which catalyses the reaction between diisopropylbenzole and 3,3',5,5' tetramethylbenzidine. This results in changes of colour ranging from orange through green to dark blue. Every green or greenish spot that appears in the reaction field within one minute is significant; the urine sample should be investigated more closely. In the case of menstruating women, blood is often, but not always, found in the urine. The significance of the detection of minor traces of blood in the urine varies from patient to patient, so that these samples need to be clinically investigated.

**Nitrite:**  
This test is based on the conversion of nitrate into nitrite owing to the presence of gram-negative bacteria in the urine. In an acidic environment the nitrite in the urine reacts with P-arsanilic acid to produce a diazonium compound. This compound then combines with 1 N-(1-naphthyl)-ethylene- diamine and produces pink coloration. Nitrite is not detectable in normal urine. With certain infections the nitrite field is positive, depending on how long the urine was held in the bladder before testing. Only 40 % of positive cases are detected by means of the nitrite test if the urine was only incubated in the bladder for a short period; however, the detection rate rises to 80 % if the urine was in the bladder for longer than 4 hours.

**Protein:**  
This reaction is based on a phenomenon known as the "protein error" of pH indicators; a highly buffered indicator changes its colour in the presence of proteins (anions) as a result of the indicator releasing hydrogen ions into the protein. When the pH value is constant, any green colour is due to the presence of protein. The coloration ranges from yellow to yellowish green in the case of negative results and from green to greenish blue in the case of positive results. 1 - 14 mg/dl protein can be excreted by healthy kidneys. Any coloration of the test field that is more intensive than the reference field is a clear indication of proteinuria. If urine has a high density, the colour on the test field may match that on the colour chart even with normal protein concentrations. Medical advice should be sought in order to evaluate the significance of traces of protein.

**LIMITATIONS OF THE PROCEDURE**

It is possible for this test to give false results. Consult your doctor before taking any medicines.

**Leukocytes:**  
The results should not be read off until 60 - 120 seconds after testing, so as to allow the colours to develop fully. The stronger the colour change, the more leukocytes are present in the urine sample. High density urine or a heavy concentration of glucose (≥ 2000 mg/dl) can lead to artificially low results. The presence of cephalixin, cephalotin or oxalic acid in high concentrations can also lead to artificially low results. Tetracyclin can lead to reduced reactivity; drugs in high dosages can lead to false negative results. When a high proportion of protein is present, the coloration may be weaker. This test does not react to erythrocytes or bacteria which are normally present in urine.

**Blood:**  
A uniform blue colour indicates the presence of myoglobin, haemoglobin or haemolysed erythrocytes. Compact or scattered blood clots indicate intact erythrocytes. For more precise test results, separate colour charts for haemoglobin and erythrocytes are included in the test kit. Positive urine test results often occur in women during menstruation. It should be noted that urine with a raised pH-value reduces sensitivity; on the other hand, a high or moderate concentration of ascorbic acid may prevent the colour change. Microbial peroxidase in association with a urinary tract infection may cause a false positive test result. The test reacts somewhat more strongly to haemoglobin and free myoglobin than to intact erythrocytes.

**Nitrite:**  
The test is nitrite specific and does not react to any substances normally excreted in urine. Any shade of uniform pink or red colour should be regarded as a positive result and indicates the presence of nitrite. The intensity of the colour is not proportional to the number of bacteria in the urine sample. Pink spots or pink-coloured edges should not be interpreted as a positive test result. If the reaction field that has changed colour is held against a white background, a relatively low nitrite content may be revealed which might otherwise go undetected. An ascorbic acid content below 30 mg/dl can lead to false negative results in urine samples which contain less than 0.05 mg/dl of nitrite ions. The sensitivity of the test is reduced in urine which is extremely alkaline. To obtain precise results, the patient must have stopped taking antibiotics at least three days before the test is performed. A negative result never rules out the possible presence of bacteria. Negative results can even occur in the case of urinary tract infections caused by other organisms which do not contain reductases that convert nitrate into nitrite, or if the urine has not remained in the bladder for a sufficient length of time (at least 4 hours) for the nitrate to be converted into nitrite, or if the patient follows a low-salt diet.

**Proteins:**  
Any green coloration indicates the presence of protein in the urine. The test shows protein in low concentrations and reacts less strongly to haemoglobin, globulin and mucu-proteins. Even if the test result is negative, other proteins may be present. False positive results can be obtained in the case of buffered urine or alkaline urine. Urine samples which are contaminated with quaternary ammonium salt or washing solutions containing chlorhexidine may cause false positive results. Urine samples with a high specific density can lead to false negative results.

**PERFORMANCE CHARACTERISTICS**

The effectiveness of the test strips for urine tests has been investigated in laboratory and clinical studies. Several factors are of decisive importance for interpreting the visual results: variations in the perception of colours, the presence or absence of retardants and the light conditions during evaluation of the test.

**BIBLIOGRAPHY**

- Tietz NW. Clinical Guide to Laboratory Tests. W.B. Saunders Company. 1976.
- Burtis CA, Ashwood ER. Tietz Textbook of Clinical Chemistry 2nd Ed. 2205, 1994.

**LIST OF SYMBOLS**

	Important: follow instructions		Number of test strips per test kit		Manufacturer
	In-vitro diagnostic agent		Expiry date		For single use only
	Min. and max. temperatures 2 ° - 30 °C		Lot number		Catalogue no.

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Geratherm Medical AG